



ergonomic products
Performance Without Compromise™

Serial Number Locations

Patient Chairs

Move the chair up and the serial number is located below the bottom seat on the left side of the seating position.

EP Chair



ETL Chair




Workstation

Serial number is located on the back of the unit next to the water bottle.



Ortho Cart

Serial number is located inside the unit under the corian top.



Hygiene Workstation

Serial number is located on the back of the unit.



clientsupport@desergo.com
866.374.6487




ergonomic products
Performance Without Compromise™

Serial Number Locations

Small Room Solutions (SRS) Workstation / Compact Cart

Serial number is located under the unit.



Spring Health Curing Light

Serial number is printed on the body (5 digit number.)



Parkell Scaler

Serial number is located on the bottom of the unit (5 digit number)



Bien Air Motor

Serial number is engraved on the body.



clientsupport@desergo.com
866.374.6487



ergonomic products
Performance Without Compromise™

Serial Number Locations

L3 TBL Track Light



Serial number is located on the metal box, behind the monitor.



L4 Aurora Track Light



Serial number is located inside the track.



clientsupport@desergo.com
866.374.6487



ergonomic products
Performance Without Compromise™

Serial Number Locations

L5 Fast Track Light

Serial number is located inside the track.



Track Light Illuminator

Serial number is on the side of the light metal cover.



clientsupport@desergo.com
866.374.6487